




Kevin Sebley - Endoscopy

 VisionScope Technologies



ARTHROSCOPY



>> Arthroscopy.

Arthroscopes		l=110 mm	l=175 mm	l=187,5 mm
Ø 2,7 mm	0°	30-0001-00	30-0004-00	30-0007-00
Ø 2,7 mm	30°	30-0002-00	30-0005-00	30-0008-00
Ø 2,7 mm	70°	30-0003-00	30-0006-00	30-0009-00
Ø 4,0 mm	0°		30-0010-00	
Ø 4,0 mm	30°		30-0011-00	
Ø 4,0 mm	45°		30-0085-00	
Ø 4,0 mm	70°		30-0012-00	

autoclavable (134°C / 273°F)

All sizes of arthroscopes can be supplied compatible with WOLF, AESCULAP
STRYKER, OLYMPUS, STORZ and DYONICS.

Kevin Sebley - Endoscopy



Web: www.arthroscopyuk.co.uk

Email: visionscope@live.co.uk

Tel: 07761 674499

