

# VST-XM-L2 LED ENDOSCOPE HANDLE



The built-in LED bulb provides the colour temperature of approx. 6500 k corresponding to sunlight. The latest circuit technology provides optimum performance and functionality. Compact and portable, the AAA batteries allow the VST-XM-L2 to work more than 2 hours. We design and manufacture it for use either with a rigid or flexible endoscope. 5 light mode settings by extending the handle body & strobe control via the on/off button.

## TECHNICAL DATA

**Battery:** 3 x AAA type battery

**Bulb:** LED bulb, approx.10 W

**Safety standard:** IEC601-1

**Color temp.:** approx.6,500 K

**Illumination:** 1200 lumen

**Life cycle of bulb:** approx.100k hours

**Storage temp.:** -15°C to +50°C

**Operating temp.:** 0°C to 45°C

**Relative humidity:**30% to 90%

**Connection:** Female ACMI/BS.

Fits all major borescopes.

