

XD-303 LED LIGHT SOURCE



Features:

This product uses 20W LED with high luminous efficiency as the light source.

Long bulb life; small product size,

Light weight, handle configuration, easy to carry.

Widely used in:uterine cavity,anorectal,otolaryngology examination

Specification:

Rated voltage: AC90~240V 50Hz-60Hz

Bulb power: LED 20W

Bulb color temperature: 5000±500K

Lamp life: 20000h

Image plane illumination: 1,600,000Lux

Net weight: 1.6 Kg

